


Example: PLT_CAR_L5_AC_A00



Type of carrier	Construction	Orientation	Number of..	Labware info	Revision
PLT Plate carrier	CAR Standard Carrier APE Application engineering DAT Deck adaptor template	L Landscape P Portrait	Plate positions: 3 4 5	AC Deepwell plate (Archive) MD Medium Density (96/384) HD High Density (1536 well) LI Limbro ST Stack	A00, A01
SMP Sample carrier	CAR Standard		Tubes: 12 16 24 32	15x75 Tube size	A00, A01
TIP Tip carrier	CAR Standard Carrier	L Landscape P Portrait	1000ul Channel: 288 384 480 5mlTips: 72 96 120 384 Head: 1920	LT Low volume 10ul 50ul 50ul Tip ST Standard vol. 300ul HT High volume 1000ul Filter 5mlT 5ml Tip Size BC = <i>Barcoded Tip rack</i> NTR = <i>Nestable Tip Rack</i>	A00, A01
RGT Reagent	CAR Standard		Reagent troughs: 3, 4, 5	R Reagent	A00, A01
CTR Control carrier				C Controls	A00, A01
VER Verification					A00, A01